

Product Datasheet

IL2RA/CD25 Mouse Monoclonal Antibody (PE-Cy5) (orb2284236)

| | |
|-----------------------------|---|
| Catalog Number | orb2284236 |
| Category | Antibodies |
| Description | IL2RA/CD25 Mouse Monoclonal Antibody (PE-Cy5) is a PE/Cy5 conjugated antibody targeting IL2RA. This antibody is suitable for FC. It exhibits reactivity with Human samples. |
| Target | IL2RA |
| Clonality | Monoclonal |
| Species/Host | Mouse |
| Isotype | IgG |
| Conjugation | PE/Cy5 |
| Reactivity | Human |
| Concentration | 1mg/ml |
| Buffer/Preservatives | 0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol. |
| Purification | Affinity purified by Protein G |
| Immunogen | KLH conjugated synthetic peptide derived from human IL2RA/CD25 |
| Tested applications | FC |
| Dilution range | Flow-Cyt=1ug/Test |
| Antibody Type | Primary Antibody |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

| | |
|------------------------|---|
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |
| Note | For research use only |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)