

## Product Datasheet

### CCR8 (Human) Recombinant Protein (P01) (orb2289713)


|                            |                                                                                                                                                                                                                                                                                                                                                                                      |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human CCR8 full-length ORF ( NP_005192.1, 1 a.a. - 355 a.a.) recombinant protein with GST-tag at N-terminal.                                                                                                                                                                                                                                                                         |
| <b>Tested Applications</b> | AP, Array, ELISA, WB                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Preservatives</b>       | 50 mM Tris-HCl, 10 mM reduced Glutathione, pH=8.0 in the elution buffer.                                                                                                                                                                                                                                                                                                             |
| <b>Storage</b>             | Store at -80°C. Aliquot to avoid repeated freezing and thawing.                                                                                                                                                                                                                                                                                                                      |
| <b>Tag</b>                 | GST                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                                                                                                                                |
| <b>Application notes</b>   | Best use within three months from the date of receipt of this protein.                                                                                                                                                                                                                                                                                                               |
| <b>Protein Sequence</b>    | MDYTLDSLVTTVTDYYPDIFSSPCDAELIQTN GKLLAVFYCLLFVFSLLGNSLVILVLVVK<br>KLRSITDVYLLNLALSDLLVFVSFPFQTYLLDQWVFGTVMCKVVSGFYIIGFYSSMFFITLMS<br>VDRLAVVHAVYALKVRTIRMGTTLC LAVWLTAIMATIPLL VFYQVASEDGV LQCYSFYNNQ<br>TLKWKIFTNFKMNLGLLPFTIFMFCYIKILHQLKRCQNHNKTKAIRLV LIVVIASLLFWVPFNV<br>VLFLTSLHSMHILDGCSISQQLTYATHVTEIISFTHCCVNPVIYAFVGEKFKKHLSEIFQKSCSQ<br>IFNYLGRQMPRESCEKSSSCQQHSSRSSVDYIL |
| <b>Clonality</b>           | Recombinant                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Antibody Type</b>       | Recombinant Antibody                                                                                                                                                                                                                                                                                                                                                                 |
| <b>NCBI</b>                | <b>NP_005192.1</b>                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Dilution Range</b>      | Heating may cause protein aggregation. Please do not heat this product before electrophoresis.                                                                                                                                                                                                                                                                                       |
| <b>Expiration Date</b>     | 6 months from date of receipt.                                                                                                                                                                                                                                                                                                                                                       |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



---

**Biorbyt Ltd.**

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

**Biorbyt LLC.**

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)