

Product Datasheet

PFN1 polyclonal antibody (A01) (orb2292420)


Description	Mouse polyclonal antibody raised against a full-length recombinant PFN1.
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	PFN1 (AAH06768, 1 a.a. ~ 140 a.a) full-length recombinant protein with GST tag.
Preservatives	50 % glycerol
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	MAGWNAYIDNLMADGTCQDAAIVGYKDSPSVWAAVPGKTFVNITPAEVGVLVGKDRSSFY VNGLTLLGGQKCSVIRDSLLQDGEFSMDLRTKSTGGAPT FNVTVTKTDKTLVLLMGKEGVH GGLINKKCYEMASHLRRSQY
Clonality	Polyclonal
NCBI	AAH06768
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



PFN1 polyclonal antibody (A01), Western Blot analysis of PFN1 expression in HeLa.



PFN1 polyclonal antibody (A01), Western Blot analysis of PFN1 expression in Raw 264.7.



Western Blot detection against Immunogen (41.51 KDa).

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)