

Product Datasheet

FHL2 monoclonal antibody (M01), clone 2G3-1A5 (orb2293037)

Catalog Number	orb2293037
Category	Antibodies
Description	Mouse monoclonal antibody raised against a full length recombinant FHL2.
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1 kappa
Conjugation	Unconjugated
Reactivity	Human
Buffer/Preservatives	In 1x PBS, pH 7.4
Immunogen	FHL2 (AAH14397, 1 a.a. ~ 279 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.
Protein Sequence	MTERFDCHHCNESLFGKKYILREESPYCVVCFETLFAANTCEECKGKPIGCDCKDLSYKDRHW HEACFHCSQCRNSLVDKPFAAKEDQLLCDYCYSNEYSSKQCCKKTIMPGRKMEYKGSS WHETCFICHRCQQPIGTKSFIPKDNQNFVPCYEKQHAMQCVQCKKPITGGVTYREQPW HKECFVCTACRKQLSGQRFTARDDFAYCLNCFCDLYAKKCAGCTNPISGLGGTKYISFEER QWHNDCFNCKKCSLSLVGRGFLTERDDILCPDCGKDI
Tested applications	ELISA, IF, WB
Clone Number	2G3-1A5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Note For research use only


NCBI **AAH14397**



FHL2 monoclonal antibody (M01), clone 2G3-1A5 Western Blot analysis of FHL2 expression in HeLa.



Immunofluorescence of monoclonal antibody to FHL2 on HeLa cell. [antibody concentration 10 ug/ml]



Western Blot detection against Immunogen (56.43 KDa).

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)