

Product Datasheet

EMR1 purified MaxPab mouse polyclonal antibody (B01P) (orb2293243)

Description	Mouse polyclonal antibody raised against a full-length human EMR1 protein.
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	EMR1 (AAH59395, 1 a.a. ~ 821 a.a) full-length human protein.
Preservatives	In 1x PBS, pH 7.4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	MRGFNLLLFWGCCVMHSEWEGHIRPTRKPNTKGNNCRDSTLCPAYATCTNTVDSYYCACK QGFLSSNGQNHFKDPGVRCKDIDECSQSPQPCGPNSSCKNLSGRYKCSCLDGFSSPTGND WVPGKPGNFSCDINECLTSSVCPEHSDCVNSMGSYSCSQVGFISRNSTCEDVDECADP RACPEHATCNNTVGNYSFCNPGFESSGHLSFQGLKASCEDIDECTEMCPINSTCTNTPG SYFCTCHPGFAPSNGQLNFTDQGVCECRDIDECRQDPSTCGPNSICTNALGSYSCGCIAGFH PNPEGSQKDGNFSCQRVLFKCKEDVIPDNKQIQQCQEGTAVKPAYVSFCAQINNIFSVLDK VCENKTTVSLKNTTESFVPVLKQISTWTKFTKEETSSLATVFLESVESMTLASFWKPSANIT PAVRTEYLDIESKVINKECSEENVTLDLVAKGDKMKIGCSTIEESESTETTGVAFVSFVGMES VLNERFFKDHQAPLTTSEIKLKMNSRVVGGIMTGEKKDGFSDPIIYTLENIQPKQKFERPICVS WSTDVKGGRWTSFGCVILEASETYTICSCNQMANLAIMASGELTMGCAIAGFLHYLFLACF FWMLVEAVILFLMVRNLKVNYFSSRNKMLHICAFGYGLPMLVVVISASVQPQGYGMHNR CWLNTETGFIWSFLGPVCTVIVINSLLLTWTLWILRQRLSSVNAEVSTLKDTRLLTFKAFAQL FILGCSWVLGIFQIGPVAGVMAYLFTIINSLQGAFIFLIHCLLNGQVREEYKRWITGKTKPSSQ SQTSRILLSSMPSASKTG
Clonality	Polyclonal

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

NCBI**AAH59395****Expiration Date**

12 months from date of receipt.



FACS analysis of negative control 293 cells (Black) and EMR1 expressing 293 cells (Green) using EMR1 purified MaxPab mouse polyclonal antibody.



Western Blot analysis of EMR1 expression in transfected 293T cell line by EMR1 MaxPab polyclonal antibody. Lane 1: EMR1 transfected lysate (90.31 KDa). Lane 2: Non-transfected lysate.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)