

Product Datasheet

CSF3 MaxPab rabbit polyclonal antibody (D01) (orb2294046)


Description	Rabbit polyclonal antibody raised against a full-length human CSF3 protein.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IP, WB
Immunogen	CSF3 (NP_000750.1, 1 a.a. ~ 207 a.a) full-length human protein.
Preservatives	No additive
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	MAGPATQSPMKLMALQLLLWHSALWTVQEATPLGPASSLPQSFLKCLEQVRKIQGDGAA LQEKLVSECATYKLCHPEELVLLGHSLGIPWAPLSSCPSQALQLAGCLSQLHSGFLYQGGL QALEGISPELGPTLDTLQLDVADFATTIWQQMEELGMAPALQPTQGAMPAFASAFQRRAG GVLVASHLQSFLEVSYRVLRLHAQP
Clonality	Polyclonal
NCBI	NP_000750.1
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Immunofluorescence of purified MaxPab antibody to CSF3 on HeLa cell. [antibody concentration 10 ug/ml]



Immunoprecipitation of CSF3 transfected lysate using anti-CSF3 MaxPab rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with CSF3 purified MaxPab mouse polyclonal antibody (B01P) (orb2294047).



Western Blot analysis of CSF3 expression in transfected 293T cell line by CSF3 MaxPab polyclonal antibody. Lane 1: CSF3 transfected lysate (22.3 KDa). Lane 2: Non-transfected lysate.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)