

Product Datasheet

CNTF purified MaxPab mouse polyclonal antibody (B01P) (orb2294283)

Catalog Number	orb2294283
Category	Antibodies
Description	Mouse polyclonal antibody raised against a full-length human CNTF protein.
Clonality	Polyclonal
Species/Host	Mouse
Conjugation	Unconjugated
Reactivity	Human
Buffer/Preservatives	In 1x PBS, pH 7.4
Immunogen	CNTF (NP_000605.1, 1 a.a. ~ 200 a.a) full-length human protein.
Protein Sequence	MAFTEHSPLTPHRRDLCRSIWLARKIRSDLTALTESYVKHQGLNKNINLDSADGMPVASTD QWSELTEAERLQENLQAYRTFHVLLARLLEDQQVHFTPTEGDFHQAIHTLLLQVAAFAYQIE ELMILLEYKIPRNEADGMPINVGDGGLFEKKLWGLKVLQELSQWTVRSIHDLRFISSHQTGIP ARGSHYIANNKKM
Tested applications	WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
NCBI	NP_000605.1

Biorbyt Ltd.


7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558



Western Blot analysis of CNTF expression in transfected 293T cell line by CNTF MaxPab polyclonal antibody. Lane 1: CNTF transfected lysate(22 KDa). Lane 2: Non-transfected lysate.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)