

Product Datasheet

CBL purified MaxPab rabbit polyclonal antibody (D01P) (orb2295022)

Description	Rabbit polyclonal antibody raised against a full-length human CBL protein.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	PLA, WB
Immunogen	CBL (NP_005179.2, 1 a.a. ~ 906 a.a) full-length human protein.
Preservatives	In 1x PBS, pH 7.4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	MAGNVKKSSGAGGGSGSGSGSGGLIGLMKDAFQPHHHHHHLSPHPPGTVDKKMVEK CWKLMDKVVRLCQNPKLALKNSPPYILDLLPDTYQHLRTILSRYEGKMETLGENEYFRVFME NLMKKTQTISLFKEGKERMYEENSQPRRNLTKLSLIFSHMLAELKGIFPSGLFQGDTFRITK ADAAEFWRKAFGEKTIVPWKSFRQALHEVHPISSGLEAMALKSTIDLTCNDYISVFEDIFTR LFQPWSSLLRNWNSLAVTHPGYMAFLTYDEVKARLQKFIHKPGSYIFRLSCTRLGQWAIGYV TADGNILQTIPH NKPLFQALIDGFREGFYLPDGRNQNPDLTGLCEPTQDHIKVTQEYELY CEMGSTFQLCKICAENDKDVKIEPCGHLMTSCLTSWQSEGGQGPCFCRCEIKGTEPIVVD PFDPRGSGSLLRQGAEGAPSPNYDDDDDERADDTLFMMKELAGAKVERPPSPFSMAPQAS LPPVPPRLDLLPQRVCVPSSASALGTASKAASGSLHKDKPLVPPTLRDLPPPPPPDRPYSVG AESRPQRRPLPCTPGDCPSRDKLPPVPSSRLGDSWLPRLPIPKVPVSAPSSSDPWTGRELTN RHSLPFSLSQMEPRPDVPRLGSTFSLDTSMSMNSSPLVGPECDHPKIKPSSSANAIYSLAA RPLVPKLPPEGQCEGEEDTEYMTPSSRPLRPLDTSQSSRACDCDQQIDSCITYEAMYNIS QAPSITESSTFGEGNLAAAHANTGPEESENEDDGYDVPKPPVPAVLARRTLSDISNASSSFG WLSLDGDPTTNVTEGSQVPERPPKPFPRRINSEKAGSCQQGSGPAASAATASPQLSSEIE NLMSQGYSYQDIQKALVIAQNNIEMAKNILREFVSISSPAHVAT


Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+4414156518558)


Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)


Clonality	Polyclonal
NCBI	NP_005179.2
Expiration Date	12 months from date of receipt.



CBL MaxPab rabbit polyclonal antibody. Western Blot analysis of CBL expression in MCF-7.



CBL MaxPab rabbit polyclonal antibody. Western Blot analysis of CBL expression in mouse spleen.




Proximity Ligation Analysis of protein-protein interactions between CBL and PLCG1. HeLa cells were stained with anti-CBL rabbit purified polyclonal 1:1200 and anti-PLCG1 mouse monoclonal antibody 1:50. Each red dot represents the detection of protein-protein interaction complex, and nuclei were counterstained with DAPI (blue).

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)



Western Blot analysis of CBL expression in transfected 293T cell line by CBL MaxPab polyclonal antibody. Lane 1: CBL transfected lysate(99.60 KDa). Lane 2: Non-transfected lysate.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)