

Product Datasheet

ATF1 purified MaxPab rabbit polyclonal antibody (D01P) (orb2295751)


Description	Rabbit polyclonal antibody raised against a full-length human ATF1 protein.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	ATF1 (NP_005162.1, 1 a.a. ~ 271 a.a) full-length human protein.
Preservatives	In 1x PBS, pH 7.4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	MEDSHKSTTSETAPQPGSAVQGAHISHIAQQVSSLSESEESQDSSDSIGSSQKAHGILARRP SYRKILKDLSSDTRGRKGDGENSGVSAAVTSMSPVTPPIYQTSSGQYIAIAPNGALQLASPG TDGVQGLQTLTMTNSGSTQQGTTILQYAQTS DGQQILVPSNQVVVQTASGDMQTYQIRTT PSATSLPQTVVMTSPVTLTSQTTKTDDPQLKREIRLMKNREAARECRRKKKEYVKCLENRVA VLENQNKT LIEELKTLKDLYSNKSV
Clonality	Polyclonal
NCBI	NP_005162.1
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



ATF1 MaxPab rabbit polyclonal antibody. Western Blot analysis of ATF1 expression in HeLa.



ATF1 MaxPab rabbit polyclonal antibody. Western Blot analysis of ATF1 expression in mouse brain.




ATF1 MaxPab rabbit polyclonal antibody. Western Blot analysis of ATF1 expression in NIH/3T3.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



ATF1 MaxPab rabbit polyclonal antibody. Western Blot analysis of ATF1 expression in PC-12.



ATF1 MaxPab rabbit polyclonal antibody. Western Blot analysis of ATF1 expression in Raw 264.7.



Western Blot analysis of ATF1 expression in transfected 293T cell line by ATF1 MaxPab polyclonal antibody. Lane 1: ATF1 transfected lysate(29.20 KDa). Lane 2: Non-transfected lysate.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)