

Product Datasheet

ANXA5 MaxPab rabbit polyclonal antibody (D01) (orb2296045)


Description	Rabbit polyclonal antibody raised against a full-length human ANXA5 protein.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	ANXA5 (NP_001145.1, 1 a.a. ~ 320 a.a) full-length human protein.
Preservatives	No additive
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Protein Sequence	MAQVLRGTVTDFPGFDERADAETLRKAMKGLGTDEESILTLTSSRNAQRQEISAAFKTLFG RDLLDDLKSELTGKFEKLIVALMKPSRLYDAYELKHALKGAGTNEKVLTEIIASRTPEELRAIK QVYEEYGGSSLEDDVVGDTSGYYQRMLVVLLQANRDPDAGIDEAQVEQDAQALFQAGELK WGTDEEKFITIFGTRSVSHLRKVFDKYMTISGFQIEETIDRETSGNLEQLLLAVVKSIRSIPAYL AETLYYAMKGAGTDDHTLIRVMVSRSEIDLFNIRKEFRKNFATSLYSMIKGDTSGDYKKALLL LCGEDD
Clonality	Polyclonal
NCBI	NP_001145.1
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.


68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



ANXA5 MaxPab rabbit polyclonal antibody. Western Blot analysis of ANXA5 expression in HepG2.



ANXA5 MaxPab rabbit polyclonal antibody. Western Blot analysis of ANXA5 expression in human liver.




ANXA5 MaxPab rabbit polyclonal antibody. Western Blot analysis of ANXA5 expression in mouse liver.

Biorbyt Ltd.


7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Immunoprecipitation of ANXA5 transfected lysate using anti-ANXA5 MaxPab rabbit polyclonal antibody and Protein A Magnetic Bead, and immunoblotted with ANXA5 purified MaxPab mouse polyclonal antibody (B01P) (orb2296046).



Western Blot analysis of ANXA5 expression in transfected 293T cell line by ANXA5 MaxPab polyclonal antibody. Lane 1: ANXA5 transfected lysate(35.9 KDa). Lane 2: Non-transfected lysate.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)