

## Product Datasheet

### CENPC Rabbit monoclonal antibody (orb2297613)

|                            |                                                                                                                                         |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog Number</b>      | orb2297613                                                                                                                              |
| <b>Category</b>            | Antibodies                                                                                                                              |
| <b>Description</b>         | CENPC Rabbit monoclonal antibody                                                                                                        |
| <b>Clonality</b>           | Monoclonal                                                                                                                              |
| <b>Species/Host</b>        | Rabbit                                                                                                                                  |
| <b>Conjugation</b>         | Unconjugated                                                                                                                            |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                       |
| <b>Form/Appearance</b>     | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.                                                                |
| <b>Purification</b>        | Affinity-chromatography                                                                                                                 |
| <b>Immunogen</b>           | A synthesized peptide derived from human CENPC                                                                                          |
| <b>UniProt ID</b>          | <b>Q03188</b>                                                                                                                           |
| <b>MW</b>                  | 107kDa                                                                                                                                  |
| <b>Tested applications</b> | FC, ICC, IHC, IP, WB                                                                                                                    |
| <b>Dilution range</b>      | WB 1:500~1:2000 IHC 1:100~1:500 ICC/IF 1:50~1:200 IP 1:50 FC 1:200                                                                      |
| <b>Antibody Type</b>       | Primary Antibody                                                                                                                        |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only                                                                                                                   |

#### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

#### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)