

Product Datasheet

Methyl ganoderenate F (orb2299288)

Catalog Number	orb2299288
Category	Small Molecules
Description	Methyl ganoderenate F is a high-purity reference standard suitable for analytical research. It is utilized in both in vitro and in vivo studies to investigate the bioactive triterpenoids found in Ganoderma lucidum, supporting pharmacological and phytochemical research.
Target	Others
MW	524.7
Formula	$C_{31}H_{40}O_7$
Storage	-20°C
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)