

Product Datasheet

Ganoderic acid GS-3 (orb2299310)

Description	Ganoderic acid GS-3 is a natural product that can be used as a reference standard. The CAS number of Ganoderic acid GS-3 is 1206781-66-9.
Storage	-20°C
Note	For research use only
Formula	$C_{32}H_{46}O_8$
Purity	98.00%
MW	558.7
CAS Number	1206781-66-9
SMILES	<chem>C[C@@]12[C@](C)(C3=C(C(=O)[C@H]1OC(C)=O)[C@]4(C)[C@@](C[C@@H]3O)(C(C)(C)[C@@H](O)CC4)[H])C(=O)C[C@@]2([C@@H](CC/C=C(/C(O)=O)\C)C)[H]</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)