

Product Datasheet

Rabbit anti-FOXL2 / BPES antibody, clone SQab30327 (monoclonal) (orb2302742)

Catalog Number	orb2302742
Category	Antibodies
Description	Rabbit anti-FOXL2 / BPES antibody, clone SQab30327 (monoclonal)
Clonality	Monoclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Purification with Protein A.
Buffer/Preservatives	PBS, 0.01% Sodium azide, 40% Glycerol and 0.05%rAlbumin.
Immunogen	Recombinant protein of Human FOXL2.
Tested applications	IHC-P
Clone Number	SQab30327
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)