

Product Datasheet

Rabbit anti-Cytokeratin (pan) antibody, clone SQab20227 (monoclonal) (orb2302792)

Catalog Number	orb2302792
Description	Rabbit anti-Cytokeratin (pan) antibody, clone SQab20227 (monoclonal)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	Synthetic peptide within aa. 100-200 of Human Type I Cytokeratins and synthetic peptide within aa. 300-400 of Human Type II Cytokeratins.
Preservatives	PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% rAlbumin.
Form/Appearance	Purification with Protein A.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Monoclonal
Clone Number	SQab20227
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)