

Product Datasheet

Rabbit Cytokeratin 17 antibody, clone SQab19157 (monoclonal) (orb2302836)

Catalog Number	orb2302836
Category	Antibodies
Description	Rabbit Cytokeratin 17 antibody, clone SQab19157 (monoclonal)
Clonality	Monoclonal
Species/Host	Rabbit
Conjugation	Unconjugated
Reactivity	Human
Buffer/Preservatives	PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% rAlbumin.
Purification	Purification with Protein A.
Immunogen	Synthetic peptide within aa. 400-500 of Human Cytokeratin 17.
Tested applications	IHC-P
Clone Number	SQab19157
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)