

## Product Datasheet

# Rabbit HA tag antibody, clone SQab18117 (monoclonal) (orb2302851)

<b>Catalog Number</b>	orb2302851
<b>Category</b>	Antibodies
<b>Description</b>	Rabbit HA tag antibody, clone SQab18117 (monoclonal)
<b>Clonality</b>	Monoclonal
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Reactivity</b>	Bovine, Canine, Gallus, Human, Monkey, Mouse, Rat
<b>Buffer/Preservatives</b>	PBS, 0.01% Sodium azide, 40% Glycerol and 0.05% rAlbumin.
<b>Purification</b>	Purification with Protein A.
<b>Immunogen</b>	YPYDVPDYA (influenza hemagglutinin-HA-epitope) conjugated to KLH.
<b>Tested applications</b>	FC, ICC, IP, WB
<b>Clone Number</b>	SQab18117
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

### Biorbyt Ltd.

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

### Biorbyt LLC

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)