

Product Datasheet

500bp DNA Ladder 500bp-5000bp (orb2303096)

Catalog Number	orb2303096
Category	Tools
Description	500 bp Ladder DNA Marker is composed of 8 DNA fragments with a length of 500 bp to 5,000 bp. dissolved in 1×Loading Buffer, and can be used for direct electrophoresis of 5-10 µl during used, which is very convenient to use. Two optimized Loading Buffer dyes are added without blocking the brightness of DNA band fluorescence.*Strips are fine-tuned, and every clip is clear enough
MW	500bp-5000bp
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)