

Product Datasheet

CY3-conjugated Goat Anti-Rabbit IgG (orb2303108)

Description	CY3-conjugated Goat Anti-Rabbit IgG
Reactivity	Rabbit
Conjugation	Cy3
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C for one year. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Absorption Peak: 550nm. Emission Peak: 570nm
Antibody Type	Secondary Antibody
Dilution Range	1:20-1:100 for most applications
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)