

Product Datasheet

Prdx1 (phospho-Tyr194) rabbit pAb Antibody (orb2304195)

Catalog Number	orb2304195
Category	Antibodies
Description	Prdx1 (phospho-Tyr194) rabbit pAb Antibody
Clonality	Polyclonal
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse, Rat
Concentration	1 mg/ml
Immunogen	Synthesized phospho peptide around human Prdx1 (Tyr194)
UniProt ID	Q06830
Tested applications	ELISA, IHC, WB
Dilution range	WB 1:500-2000; IHC-p 1:50-300; ELISA 2000-20000
Source	Rabbit
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)