

Product Datasheet

EPHA1 (Phospho-Tyr605) rabbit pAb Antibody (orb2305409)

Catalog Number	orb2305409
Category	Antibodies
Description	EPHA1 (Phospho-Tyr605) rabbit pAb Antibody
Clonality	Polyclonal
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Mouse
Concentration	1 mg/ml
Immunogen	Synthesized peptide derived from human EPHA1 (Phospho-Tyr605)
UniProt ID	P21709
Tested applications	IF, IHC, WB
Dilution range	IHC-p 1:50-200, WB 1:500-2000
Source	Rabbit
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558