

Product Datasheet

GGT1 (heavy chain, Cleaved-Gly380) rabbit pAb (orb2307419)

Catalog Number	orb2307419
Description	GGT1 (heavy chain, Cleaved-Gly380) rabbit pAb
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Immunogen	Synthesized peptide derived from human GGT1 (heavy chain, Cleaved-Gly380)
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Source	Rabbit
Dilution Range	WB 1:1000-2000 ELISA 1:5000-20000
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)