


## Product Datasheet

# Positive Control Lysate for anti- DC-SIGN / CD209 antibody (orb2308880)

<b>Catalog Number</b>	orb2308880
<b>Category</b>	Tools
<b>Description</b>	Positive Control Lysate for anti- DC-SIGN / CD209 antibody
<b>Tested applications</b>	WB
<b>Application notes</b>	Western Blot: 2µg/ml staining of Primary antibody in lysates orb2308880, Approx.50kDa band observed (calculated MW of 45.8 kDa according to NP_066978.1). Primary incubation 1 hour at room temperature. Secondary antibody at 1:40K
<b>Storage</b>	Store lysate at -20°C upon arrival. Avoid repeated freeze /thawed cycles.
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.



Western Blot: 2 µg/ml staining of Primary antibody in lysates, approx. 50kDa band observed (calculated MW of 45.8 kDa according to NP\_066978.1). Primary incubation 1 hour at room temperature. Secondary antibody is at 1:40K.

**Biorbyt Ltd.**

7 Signet Court, Swann Road  
Cambridge  
CB5 8LA  
United Kingdom

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

**Biorbyt LLC**

68 TW Alexander Drive  
Research Triangle Park  
Durham  
NC 27713  
United States

Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)