

Product Datasheet

ILDR2 Biosimilar Antibody (orb2321726)

Catalog Number	orb2321726
Category	Antibodies
Description	Anti-ILDR2(bapotulimab biosimilar) mAb
Target	ILDR2
Clonality	Monoclonal
Species/Host	Human
Isotype	IgG2(S257A)
Conjugation	Unconjugated
Reactivity	Human
Purification	Purified from cell culture supernatant by affinity chromatography
UniProt ID	Q71H61
Tested applications	FC
Dilution range	Flow Cyt 1:100
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % – 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Antibody Type	Biosimilar Antibody
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Note

For research use only



Flow cytometry analysis with 1 $\mu\text{g/mL}$ Anti-ILDR2 (bapotulimab biosimilar) mAb (orb2321726) on Expi293 cells transfected with Human ILDR2 protein (Blue histogram) or Expi293 transfected with irrelevant protein (Red histogram).

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)