

Product Datasheet

Biotinylated Anti-EPCAM(oportuzumab biosimilar) mAb (orb2321785)

Catalog Number	orb2321785
Description	Biotinylated Anti-EPCAM(oportuzumab biosimilar) mAb
Species/Host	Human Humanized
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	EPCAM
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions.
Isotype	IgG1
Clonality	Monoclonal
Antibody Type	Biosimilar Antibody
Uniprot ID	P16422
Dilution Range	ELISA 1:5000-10000, Flow Cyt 1:100
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)