

Product Datasheet

Biotinylated Anti-B7-H5(onvatilimab biosimilar) mAb (orb2321853)

Catalog Number	orb2321853
Description	Biotinylated Anti-B7-H5(onvatilimab biosimilar) mAb
Species/Host	Human
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Target	B7-H5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	IgG1
Clonality	Monoclonal
Antibody Type	Biosimilar Antibody
Uniprot ID	Q9H7M9
Dilution Range	ELISA 1:5000-10000; Flow Cyt 1:100
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)