

Product Datasheet

Biotinylated Anti-CD20 (obinutuzumab biosimilar) mAb (orb2322414)

Catalog Number	orb2322414
Description	Biotinylated Anti-CD20 (obinutuzumab biosimilar) mAb
Species/Host	Human Humanized
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Target	CD20
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8% trehalose is added as protectants before lyophilization. Please see Certificate of Analysis for specific instructions of reconstitution.
Isotype	IgG1(E239D, M241L)
Clonality	Monoclonal
Antibody Type	Biosimilar Antibody
Uniprot ID	P11836
Dilution Range	Flow Cyt 1:100
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)