

Product Datasheet

Phospho-Tau (Ser202/Thr205) Recombinant Rabbit mAb (orb2327003)

Description	Phospho-Tau (Ser202/Thr205) Recombinant Rabbit mAb
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	A synthetic phosphopeptide corresponding to residues surrounding Ser202/Thr205 of human Tau
Form/Appearance	Store at -20°C. Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% rAlbumin. Stable for 12 months from date of receipt.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Recombinant
Antibody Type	Primary Antibody
MW	50-80 kDa
Uniprot ID	P10636
Dilution Range	WB: 1:1000IP: 1:20
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)