

## Product Datasheet

### Rabbit Anti-Rat IgG H&L, BF488 conjugated (orb2328093)

|                             |                                                                         |
|-----------------------------|-------------------------------------------------------------------------|
| <b>Catalog Number</b>       | orb2328093                                                              |
| <b>Category</b>             | Antibodies                                                              |
| <b>Description</b>          | Rabbit Anti-Rat IgG H&L, BF488 conjugated                               |
| <b>Clonality</b>            | Polyclonal                                                              |
| <b>Species/Host</b>         | Rabbit                                                                  |
| <b>Isotype</b>              | IgG                                                                     |
| <b>Conjugation</b>          | BF488                                                                   |
| <b>Reactivity</b>           | Rat                                                                     |
| <b>Concentration</b>        | 2.0 mg/ml                                                               |
| <b>Buffer/Preservatives</b> | 10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol. |
| <b>Purification</b>         | Affinity purified by Protein A                                          |
| <b>Immunogen</b>            | Native Rat IgG                                                          |
| <b>Tested applications</b>  | FC, ICC, IF                                                             |
| <b>Dilution range</b>       | IF=1:100-1000, Flow-Cyt=1:100-1000, ICC/IF=1:100-1000                   |
| <b>Antibody Type</b>        | Secondary Antibody                                                      |
| <b>Storage</b>              | Store at -20°C for one year. Avoid repeated freeze/thaw cycles.         |
| <b>Note</b>                 | For research use only                                                   |

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)