

Product Datasheet

Rabbit Anti-Goat IgM, BF488 conjugated (orb2328157)

Catalog Number	orb2328157
Category	Antibodies
Description	Rabbit Anti-Goat IgM, BF488 conjugated
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	BF488
Reactivity	Goat
Concentration	2.0 mg/ml
Buffer/Preservatives	10 mM TBS (pH=7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% glycerol.
Purification	Affinity purified by Protein A
Immunogen	Native Goat IgM
Tested applications	FC, ICC, IF
Dilution range	IF=1:100-1000, Flow-Cyt=1:100-1000, ICC/IF=1:100-1000
Antibody Type	Secondary Antibody
Storage	Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)