

Product Datasheet

Recombinant human DKK-1 Protein,C-His (Active,HEK293) (orb2328311)

Catalog Number	orb2328311
Category	Proteins
Description	Recombinant human DKK-1 Protein,C-His (Active,HEK293), expressed in HEK293 cell, with a molecular weight of Approximately 35-50 kDa due to the glycosylation, and a purity of > 95% as determined by SDS-PAGE. It was purified by AC, and has a C-His tag.
Target	Dickkopf-related protein 1
Tag	C-His
Form/Appearance	Lyophilized
Buffer/Preservatives	Lyophilized after extensive dialysis against PBS.
Reconstitution	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.For long term storage it is recommended to add a carrier protein (0.1% BSA, 5% HSA, 10% FBS or 5% Trehalose).
Purity	> 95% as determined by SDS-PAGE
Purification	AC
Protein Sequence	32-265/265
UniProt ID	O94907
MW	Approximately 35-50 kDa due to the glycosylation
Expression System	HEK293 cell

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Endotoxins	1.0 EU/ μ g as determined by LAL
Biological Origin	Human
Storage	Stored at -20°C for 2 years. After reconstitution, it is stable at 4°C for 1 week or -20°C for longer (with carrier protein). It is recommended to freeze aliquots at -20°C or -70°C for extended storage.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)