

Product Datasheet

Recombinant human CCL15 Protein,N-MBP & C-His (HEK293) (orb2328349)

Catalog Number	orb2328349
Description	Recombinant human CCL15 Protein,N-MBP & C-His (HEK293)
Target	C-C motif chemokine ligand 15
Preservatives	150 mM PBS
Form/Appearance	Liquid
Storage	Stored at -70°C or -20°C for one year. Avoid repeated freeze/thaw cycles.Stored at 2 to 8 °C for 2 days after thawing.
Tag	N-MBP & C-His
Note	For research use only
Protein Sequence	22-113/113
Purity	> 85% as determined by SDS-PAGE
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)