

Product Datasheet

Recombinant human LTA Protein,C-His (HEK293) (orb2328549)

Catalog Number	orb2328549
Category	Proteins
Description	Recombinant human LTA Protein,C-His (HEK293), expressed in HEK293 cell, and a purity of > 85% as determined by SDS-PAGE. It was purified by AC, and has a C-His tag.
Target	lymphotoxin alpha
Tag	C-His
Form/Appearance	Liquid
Buffer/Preservatives	150 mM PBS
Purity	> 85% as determined by SDS-PAGE
Purification	AC
Protein Sequence	71-200/205
Expression System	HEK293 cell
Biological Origin	Human
Storage	Stored at -70°C or -20°C for one year. Avoid repeated freeze/thaw cycles.Stored at 2 to 8 °C for 2 days after thawing.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)