

Product Datasheet

Recombinant human KITLG Protein,C-His (HEK293) (orb2328622)

Catalog Number	orb2328622
Category	Proteins
Description	Recombinant human KITLG Protein,C-His (HEK293), expressed in HEK293 cell, and a purity of > 85% as determined by SDS-PAGE. It was purified by AC, and has a C-His tag.
Target	KIT ligand
Tag	C-His
Form/Appearance	Liquid
Buffer/Preservatives	150 mM PBS
Purity	> 85% as determined by SDS-PAGE
Purification	AC
Protein Sequence	121-210/273
Expression System	HEK293 cell
Biological Origin	Human
Storage	Stored at -70°C or -20°C for one year. Avoid repeated freeze/thaw cycles.Stored at 2 to 8 °C for 2 days after thawing.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558