

Product Datasheet

Recombinant human IL5 Protein,N-hFc & C-His (HEK293) (orb2328657)

Catalog Number	orb2328657
Category	Proteins
Description	Recombinant human IL5 Protein,N-hFc & C-His (HEK293), expressed in HEK293 cell, and a purity of > 85% as determined by SDS-PAGE. It was purified by AC, and has a N-hFc & C-His tag.
Target	interleukin 5
Tag	N-hFc & C-His
Purity	> 85% as determined by SDS-PAGE
Purification	AC
Protein Sequence	20-134/134
Expression System	HEK293 cell
Storage	Stored at -70°C or -20°C for one year. Avoid repeated freeze/thaw cycles.Stored at 2 to 8 °C for 2 days after thawing.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)