

Product Datasheet

Recombinant human ACE2 protein, His & Avi (HEK293), Biotin conjugated (orb2328788)

Catalog Number	orb2328788
Category	Proteins
Description	Recombinant human ACE2 protein, His & Avi (HEK293), Biotin conjugated, with a concentration of >0.5 mg/ml, expressed in HEK293 cell, with a molecular weight of 100-110 kDa, and a purity of >95% as determined by SDS-PAGE. It was purified by AC, and has a His & Avi tag.
Target	Spike glycoprotein
Tag	His & Avi
Form/Appearance	Lyophilized or Liquid
Buffer/Preservatives	Lyophilized from 0.22um filtered solution in 2X PBS (pH7.4). Normally 5% trehalose is added as protectant before lyophilization.
Reconstitution	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 µg/ml is recommended. Dissolve the lyophilized protein in distilled water.
Purity	>95% as determined by SDS-PAGE
Purification	AC
Protein Sequence	18-740/805
MW	100-110 kDa
Expression System	HEK293 cell
Endotoxins	Not analyzed

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)