

Product Datasheet

FLAG Tag Antibody (orb234224)

Catalog Number	orb234224
Category	Antibodies
Description	FLAG Tag Antibody
Target	FLAG Tag
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG2b
Conjugation	Unconjugated
Form/Appearance	Lyophilized protein G purified in PBS pH 7.4
Buffer/Preservatives	None
Immunogen	The antibody was raised against DYKDDDDK (FLAG) synthetic peptide conjugated to KLH.
Tested applications	DOT, ELISA, IHC, IP, WB
Dilution range	If reconstituted with deionized water in 100 µl: WB 1:1000-3,000, IHC 1:500-2,000. Optimal dilution has to be determined by the user.
Specificity	The antibody recognizes the N-terminal, C-terminal or internal DYKDDDDK-tagged fusion proteins.
Clone Number	FG4R
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)