

Product Datasheet

NF kB P65 (Acetyl-Lys314/Lys315) Mouse Monoclonal Antibody (APC-Cy5.5) (orb2367736)

Catalog Number	orb2367736
Description	NF kB P65 (Acetyl-Lys314/Lys315) Mouse Monoclonal Antibody (APC-Cy5.5)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	APC/Cy5.5
Tested Applications	IF
Immunogen	KLH conjugated Synthesised acetylpeptide derived from human NFkBp65 around the acetylation site of Lys314/315 IM(Ac-K)(Ac-K)SP
Target	RELA
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Clonality	Monoclonal
Clone Number	4B11
Antibody Type	Primary Antibody
Dilution Range	IF=1:100-500
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)