

Product Datasheet

P53 (Mono Methyl Lys370) Mouse Monoclonal Antibody (APC-Cy7) (orb2368541)

Description	P53 (Mono Methyl Lys370) Mouse Monoclonal Antibody (APC-Cy7)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	APC/Cy7
Tested Applications	IF
Immunogen	KLH conjugated synthesised methylpeptide derived from human P53 around the methylation site of Mono Methyl Lys370 HL(Mono Methyl K)SK
Target	TP53
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
Clone Number	5B10
Antibody Type	Primary Antibody
Dilution Range	IF=1:100-500
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)