

Product Datasheet

Histone H3 (di methyl K36) Mouse Monoclonal Antibody (PerCP-Cy5.5) (orb2369105)

Catalog Number	orb2369105
Category	Antibodies
Description	Histone H3 (di methyl K36) Mouse Monoclonal Antibody (PerCP-Cy5.5)
Target	H3C1
Clonality	Monoclonal
Species/Host	Mouse
Isotype	IgG1
Conjugation	PerCP/Cy5.5
Reactivity	Human, Mouse, Rat
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein G
Immunogen	KLH conjugated synthesised methylpeptide derived from human Histone H3 around the methylation site of di methyl K36 GV(di Methyl-K)KP
Tested applications	IF
Dilution range	IF=1:100-500
Antibody Type	Primary Antibody
Clone Number	3F7

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Storage

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)