

Product Datasheet

Histone H3 (tri methyl K79) Mouse Monoclonal Antibody (APC-Cy5.5) (orb2369300)

Description	Histone H3 (tri methyl K79) Mouse Monoclonal Antibody (APC-Cy5.5)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	APC/Cy5.5
Tested Applications	IF
Immunogen	KLH conjugated synthesised methylpeptide derived from human Histone H3 around the methylation site of Tri Methyl K79 DF(tri methyl K)TD
Target	H3C1
Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Clone Number	3C7
Antibody Type	Primary Antibody
Dilution Range	IF=1:100-500
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)