

Product Datasheet

FKBP1A Antibody (Biotin) (orb239655)

Catalog Number	orb239655
Category	Antibodies
Description	This FKBP1A Antibody, Biotin conjugated is an Biotin polyclonal product. It targets FKBP1A using a Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP1A protein (2-103AA) as the immunogen. This antibody is suitable for ELISA. Purification: >95%, Protein G purified.
Target	FKBP1A
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Biotin
Reactivity	Human
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	>95%, Protein G purified
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP1A protein (2-103AA)
UniProt ID	P62942
Tested applications	ELISA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)


Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Note For research use only

Expiration Date 12 months from date of receipt.



Western blot analysis of (Lane 1): EC109 whole cell lysate(Lane 2): 293T whole cell lysateusing Peptidyl-prolyl cis-trans isomerase FKBP1A antibody (Biotin)

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)