

Product Datasheet


AOX1 Antibody (orb239836)

Catalog Number	orb239836
Category	Antibodies
Description	This AOX1 Antibody is an unconjugated polyclonal product. It targets Sporamin A using a Recombinant Human Aldehyde oxidase protein (236-421AA) as the immunogen. This antibody is suitable for ELISA, IF, IHC. Purification: >95%, Protein G purified.
Target	Sporamin A
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	>95%, Protein G purified
Immunogen	Recombinant Human Aldehyde oxidase protein (236-421AA)
UniProt ID	Q06278
Tested applications	ELISA, IF, IHC
Dilution range	IHC-P:1:20-200

Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Biorbyt Ltd.	Biorbyt LLC
7 Signet Court, Swann Road	68 TW Alexander Drive
Cambridge	Research Triangle Park
CB5 8LA	Durham
United Kingdom	NC 27713
Email: info@biorbyt.com , support@biorbyt.com	United States
Phone: +44 (0)1223 859353 Fax: +1 (415) 651-8558	Email: info@biorbyt.com , support@biorbyt.com
	Phone: +1 (415) 906-5211 Fax: +1 (415) 651-8558

Note

For research use only



Immunohistochemical staining of human colon cancer using
Aldehyde oxidase antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)