

Product Datasheet


LDHA Antibody (orb239997)

Catalog Number	orb239997
Category	Antibodies
Description	This LDHA Antibody is an unconjugated polyclonal product. It targets Hemoglobin using a Recombinant Human L-lactate dehydrogenase A chain protein (5-323AA) as the immunogen. This antibody is suitable for ELISA, IF, IHC, WB. Purification: >95%, Protein G purified.
Target	Hemoglobin
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	>95%, Protein G purified
Immunogen	Recombinant Human L-lactate dehydrogenase A chain protein (5-323AA)
UniProt ID	P00338
Tested applications	ELISA, IF, IHC, WB
Dilution range	WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200


Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Biorbyt Ltd.	Biorbyt LLC
7 Signet Court, Swann Road	68 TW Alexander Drive
Cambridge	Research Triangle Park
CB5 8LA	Durham
United Kingdom	NC 27713
Email: info@biorbyt.com , support@biorbyt.com	United States
Phone: +44 (0)1223 859353 Fax: +1 (415) 651-8558	Email: info@biorbyt.com , support@biorbyt.com
	Phone: +1 (415) 906-5211 Fax: +1 (415) 651-8558

Note


For research use only



Immunohistochemistry of paraffin-embedded human endometrial cancer tissue using L-lactate dehydrogenase A chain antibody



Western blot aanalysis of EC109 whole cell lysate (Lane1), 293T whole cell lysate (Lane2)using L-lactate dehydrogenase A chain antibody



Western blot aanalysis of HepG2 cell using L-lactate dehydrogenase A chain antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com


Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC


68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)



Western blot analysis of HeLa whole cell lysate, A375 whole cell lysate, Jurkat whole cell lysate, 293 whole cell lysate using L-lactate dehydrogenase A chain antibody



Immunofluorescent analysis of HepG2 cells using L-lactate dehydrogenase A chain antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)