

Product Datasheet

Ca1 Antibody (Biotin) (orb241203)

Catalog Number	orb241203
Category	Antibodies
Description	This Ca1 Antibody, Biotin conjugated is an Biotin polyclonal product. It tagets Ca1 using a Recombinant Rat Carbonic anhydrase 1 protein (2-261AA) as the immunogen. This antibody is suitable for ELISA. Purification: >95%, Protein G purified.
Target	Ca1
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Biotin
Reactivity	Rat
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	>95%, Protein G purified
Immunogen	Recombinant Rat Carbonic anhydrase 1 protein (2-261AA)
UniProt ID	BOBNN3
Tested applications	ELISA
Dilution range	ELISA

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note

For research use only



Western blot analysis of (Lane 1):mouse kidney tissue (Lane 2):mouse lung tissue using Carbonic anhydrase 1 antibody (Biotin)

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)