

Product Datasheet

Phospho-AQP2 (Ser264) Rabbit Polyclonal Antibody (APC-Cy5.5) (orb2412890)

Catalog Number	orb2412890
Category	Antibodies
Description	Phospho-AQP2 (Ser264) Rabbit Polyclonal Antibody (APC-Cy5.5) is a APC/Cy5.5 conjugated antibody targeting AQP2. This antibody is suitable for ICC. It exhibits reactivity with Human samples.
Target	AQP2
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	APC/Cy5.5
Reactivity	Human
Predicted Reactivity	Bovine, Canine, Gallus, Human, Mouse, Porcine, Rat, Sheep
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human AQP2 around the phosphorylation site of Ser264 PQ(p-S)LP
Tested applications	ICC

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Dilution range	ICC/IF=1:50-200
Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)