

Product Datasheet

Phospho-JNK1 + JNK2 + JNK3 (T183 + T183 + T221) Rabbit Polyclonal Antibody (PerCP-Cy5.5) (orb2422113)

Catalog Number	orb2422113
Category	Antibodies
Description	Phospho-JNK1 + JNK2 + JNK3 (T183 + T183 + T221) Rabbit Polyclonal Antibody (PerCP-Cy5.5) is a PerCP/Cy5.5 conjugated antibody targeting MAPK8. This antibody is suitable for IF. It exhibits reactivity with Human, Mouse, Rat samples.
Target	MAPK8
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	PerCP/Cy5.5
Reactivity	Human, Mouse, Rat
Predicted Reactivity	Bovine, Canine, Porcine
Concentration	1mg/ml
Buffer/Preservatives	0.01M TBS (pH7.4) with 1% rAlbumin, 0.02% Proclin300 and 50% Glycerol.
Purification	Affinity purified by Protein A
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human JNK1/JNK2/JNK3 around the phosphorylation site of T183/T183/T221 MM(p-T)PY

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Tested applications	IF
Dilution range	IF=1:100-500
Antibody Type	Primary Antibody
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713-2847
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)