

Product Datasheet

EIF3I Antibody (orb242786)

Catalog Number	orb242786
Category	Antibodies
Description	This EIF3I Antibody is an unconjugated polyclonal product. It targets EIF3I using a Recombinant Human Eukaryotic translation initiation factor 3 subunit I protein (1-323AA) as the immunogen. This antibody is suitable for ELISA, IF, WB. Purification: >95%, Protein G purified.
Target	EIF3I
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Human, Rat
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	>95%, Protein G purified
Immunogen	Recombinant Human Eukaryotic translation initiation factor 3 subunit I protein (1-323AA)
UniProt ID	Q13347
Tested applications	ELISA, IF, WB
Dilution range	WB:1:2000-1:10000, IF:1:200-1:500

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage


Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only



Western blot analysis of EC109(Lane 1), 293T(Lane 2) whole cell lysate using Eukaryotic translation initiation factor 3 subunit I antibody



Immunofluorescence analysis of HeLa cells using Eukaryotic translation initiation factor 3 subunit I antibody

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558