

Product Datasheet

Alb Antibody (Biotin) (orb243539)

| | |
|-----------------------------|--|
| Catalog Number | orb243539 |
| Category | Antibodies |
| Description | This Alb Antibody, Biotin conjugated is an Biotin polyclonal product. It targets IgG using a Serum Albumin purified from Mouse serum (Native Protein) as the immunogen. This antibody is suitable for ELISA. Purification: >95%, Protein G purified. |
| Target | Alb |
| Clonality | Polyclonal |
| Species/Host | Rabbit |
| Isotype | IgG |
| Conjugation | Biotin |
| Reactivity | Mouse |
| Form/Appearance | Liquid |
| Buffer/Preservatives | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol |
| Purification | >95%, Protein G purified |
| Immunogen | Serum Albumin purified from Mouse serum (Native Protein) |
| UniProt ID | P07724 |
| Tested applications | ELISA |
| Dilution range | WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200 |

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Storage

Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.

Note

For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)