

Product Datasheet

RuBisCO activase Antibody (orb243651)

Catalog Number	orb243651
Category	Antibodies
Description	This RuBisCO activase Antibody is an unconjugated polyclonal product. It targets RuBisCO activase using a Ribulose-1,5-bisphosphate carboxylase oxygenase (Native Protein) as the immunogen. This antibody is suitable for ELISA, WB. Purification: >95%, Protein G purified.
Target	RuBisCO activase
Clonality	Polyclonal
Species/Host	Rabbit
Isotype	IgG
Conjugation	Unconjugated
Reactivity	Plant
Form/Appearance	Liquid
Buffer/Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Purification	>95%, Protein G purified
Immunogen	Ribulose-1,5-bisphosphate carboxylase oxygenase (Native Protein)
UniProt ID	P10871
Tested applications	ELISA, WB
Dilution range	WB:1:1000-5000

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States


Email: info@biorbyt.com, support@biorbyt.com
Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558

Storage

Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note

For research use only



Western blot analysis of Ribulose-1,5-bisphosphate carboxylase oxygenase

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: +44 (0)1223 859353 | Fax: +1 (415) 651-8558

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: +1 (415) 906-5211 | Fax: +1 (415) 651-8558